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Optilene® Mesh LP

Produktbeschreibung

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Optilene® Mesh LP

Description of the device

Optilene® Mesh LP is a mesh implant for re-inforcement of connective tissue structures. It is constructed from monofilament polypropylene, that has been knitted to a thin and elastic shape-stable mesh. After implantation the Optilene® Mesh LP adapts to the longitudinal and latitudinal expansions taking place in the connective tissue.

Optilene® Mesh LP is available undyed and with blue stripes coloured with copper phthalocyanine (Phthalocyaninato (2-) copper).

Optilene® Mesh LP does not possess any independent pharmacological properties. The polypropylene mesh is biostable and is not degraded in the body.

Indications for Use

Re-inforcement material in connective tissue weakness:

- for prosthetic hernioplasty
- for reconstruction of the chest wall
- for re-inforcement of fascial tissue

if a non-absorbable reinforcement material is required.

Optilene® Mesh LP is suitable for both conventional and laparoscopic application.

Application

Select Optilene® Mesh LP from the product range according to the size of the defect and cut it to the size to fit the defect.

When repairing inguinal hernias, the mesh should, in order to prevent recurrence, be large enough to project beyond the pubic tubercle and pass reliably round the spermatic cord at the inner ring.

Non-absorbable suture material or clips are recommended for fixing Optilene® Mesh LP.

They should be placed ≥ 10 mm from the edge of the mesh. Atraumatic round-bodied needles make suturing possible without greater damage to the implant.

Mode of action

The healing-in process is primarily influenced by the porosity of the implant. When Optilene® Mesh LP is implanted, connective tissue cells soon begin to migrate into the pores of the implant. Here there can be a brief, mild tissue reaction, which, however, supports a rapid, problem-free healing-in process. The defined pore size brings about even, tension-free penetration by the tissue. The Optilene® Mesh LP is characterized by its thin but extremely shape-stable mesh wall. It adapts optimally to the longitudinal and latitudinal motion of the patient as a result of body movement and muscle contraction.

Contraindications

Optilene® Mesh LP should not be implanted:

- in contaminated and infected areas
- in children during the growth phase

Warnings

- Direct contact between Optilene® Mesh LP and the viscera (intestines etc) must be avoided to prevent adhesion.
- Only use non-absorbable material (suture material, staples, tackers devices) together with Optilene® Mesh LP.
- Do not reuse the product: infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.

Precautions

- Before implanting the mesh the surgeon must be familiar with the surgical technique and the specific characteristics of the different polypropylene meshes.

- Optilene® Mesh LP may only be used if the packaging is undamaged and, hence, the sterility is guaranteed.

- Optilene® Mesh LP may not be used after the expiry date given on the packaging.

- Resterilization of Optilene® Mesh LP is not allowed.

Adverse events

Like all non-absorbable meshes, implantation of Optilene® Mesh LP can be associated with a limited period of local irritation in the area of the wound; a transient foreign body reaction can occasionally take the form of an inflammatory reaction. Infections and adhesions on contact with the viscera arise as a result of the intervention.

Sterilization

Optilene® Mesh LP is sterilized using ethylene oxide.

Storage

No special storage conditions are required.

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Title: FOLLETO OPTILENE MESH LP V.2 Initiator: Susana ? Calvillo

This document is signed electronically in compliance with the B. Braun electronic signature policies and procedures by following persons:

UserName: Calvillo, Susana (calvsues)
Title: SUT-PR-ES04-Ayudante Técnico Sistema Impresión
Date: Tuesday, 02 February 2016, 11:29 W. Europe Daylight Time
Meaning: Document signed as Author
=====

UserName: Garcia, Guillermo (garcgues)
Title: SU-PR-ES04-Technical support manager
Date: Wednesday, 03 February 2016, 14:40 W. Europe Daylight Time
Meaning: Precheck of Document
=====

UserName: Benade, Miguel Angel (benamies)
Title:
Date: Wednesday, 03 February 2016, 16:53 W. Europe Daylight Time
Meaning: Approve Document
=====

UserName: Piera Eroles, Eva (piermaes)
Title:
Date: Thursday, 04 February 2016, 09:21 W. Europe Daylight Time
Meaning: Approve Document
=====

UserName: Turon, Pau (turopaes)
Title: R&D, Regulatory Affairs and Quality Management Director CoE CT
Date: Thursday, 04 February 2016, 20:37 W. Europe Daylight Time
Meaning: Final Release of the Document
=====