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Contraindikations

Optilene® Mesh should not be implanted:

- in contaminated and infected areas
- in children during the growth phase

Warnings

- Direct contact between Optilene® Mesh and the viscera (intestines etc.) must be avoided to prevent adhesion.

- Only use non-absorbable material (suture material, staples, tacker devices) together with Optilene® Mesh.

- Do not reuse the product: Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.

- No special storage conditions are required.

Notes/reparative measures

- Before implanting Optilene® Mesh the surgeon must be familiar with the surgical technique and the specific characteristics of the different polypropylene meshes.

- Optilene® Mesh may only be used if the packaging is undamaged and, hence, the sterility is guaranteed.

- Optilene® Mesh may not be used after the expiry date given on the packing.

- Resterilization of Optilene® Mesh is not recommended.

Side effects

Like all non-absorbable meshes, implantation of Optilene® Mesh can be associated with a limited period of local irritation in the area of the wound: a transient foreign body reaction can occasionally take the form of an inflammatory reaction. Infections and adhesions with the viscera can occur as a result of the intervention.

Sterilization

Optilene® Mesh is sterilized using ethylene oxide.
Date of information: 04/2016

tissue cells soon begin to migrate into the pores of the implant. Here there can be a brief, mild tissue reaction, which, however, supports a rapid, problem-free healing-in process. The defined pore size brings about even, tension-free penetration by the tissue. The Optilene® Mesh is characterized by its thin but extremely shape-stable mesh wall. It adapts optimally to the longitudinal and latitudinal motion of the patient as a result of body movement and muscle contraction.

Contraindications

Optilene® Mesh should not be implanted:

- in contaminated and infected areas
- in children during the growth phase

Warnings

- Direct contact between Optilene® Mesh and the viscera (intestines etc.) must be avoided to prevent adhesion.

- Only use non-absorbable material (suture material, staples, tacker devices) together with Optilene® Mesh.

- Do not reuse the product: Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.

- No special storage conditions are required.

Notes/reparative measures

- Before implanting Optilene® Mesh the surgeon must be familiar with the surgical technique and the specific characteristics of the different polypropylene meshes.

- Optilene® Mesh may only be used if the packaging is undamaged and, hence, the sterility is guaranteed.

- Optilene® Mesh may not be used after the expiry date given on the packing.

- Resterilization of Optilene® Mesh is not recommended.

Side effects

Like all non-absorbable meshes, implantation of Optilene® Mesh can be associated with a limited period of local irritation in the area of the wound: a transient foreign body reaction can occasionally take the form of an inflammatory reaction. Infections and adhesions with the viscera can occur as a result of the intervention.

Sterilization

Optilene® Mesh is sterilized using ethylene oxide.

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(Approval & changes in printed materials for sutures)

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
**Document Control
& Signature Page**

Title: FOLLETO OPTILENE MESH V.2 Initiator: Susana ? Calvillo

This document is signed electronically in compliance with the B. Braun electronic signature policies and procedures by following persons:

UserName: Calvillo, Susana (calvsues)

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Date: Monday, 27 June 2016, 10:08 W. Europe Daylight Time

Meaning: Document signed as Author

=====

UserName: Garcia, Guillermo (garcgues)

Title: SU-PR-ES04-Technical support manager

Date: Monday, 27 June 2016, 10:33 W. Europe Daylight Time

Meaning: Precheck of Document

=====

UserName: Piera Eroles, Eva (piermaes)

Title: Senior Regulatory Affairs

Date: Thursday, 30 June 2016, 11:33 W. Europe Daylight Time

Meaning: Approve Document

=====

UserName: Turon, Pau (turopaes)

Title: R&D, Regulatory Affairs and Quality Management Director CoE CT

Date: Thursday, 30 June 2016, 16:52 W. Europe Daylight Time

Meaning: Final Release of the Document

=====