

PRODUCT DATA SHEET

Flexima® Active Drainable Bags

One-piece stoma bags for people with ileostomy or colostomy having pasty or liquid stools



The Flexima® Active range offers a **B. Braun patented skin protector** which makes you feel the difference: Its high cohesion network provides **immediate adhesivity**, **protects the skin** and presents **low risk of residue**.* The low profile skin protector of this 1-piece range is very comfortable to wear and adapts easily to all body contours giving users a **sensation of second skin****.

Flexima® Active drainable bags with integrated Roll'up clamp are available in three sizes, mini, midi and maxi.

Features

Intuitive and secure closing system:

- Outlet with integrated self-gripping soft clamp
- 7 cm wide outlet opening
- The edges are softened to increase the wearing comfort once the Roll'up outlet is closed
- Once folded and closed, the drainable bag has a shape as compact as of a colostomy bag to prevent bending effect.

High performance protected filter (HP):

- Improved deodorization capacity through its high-tech breathable membrane
- Leak-proof: water and oil resistant, permeable to gases

Colors and capacity:

In addition to the beige and transparent bags we are offering beige bags with inspection window that combine both the discretion of a beige stoma bags and the visual control to check easily stoma or stool aspect and quantity.

- Mini:
 - Bag capacity of +/- 260 ml (ISO 8670-2)
 - Cut-to-fit Ø15-45 mm: beige, transparent and beige with inspection window
- Midi:
 - Bag capacity of +/- 580 ml (ISO 8670-2)
 - Cut-to-fit Ø15-50 mm: beige, transparent and beige with inspection window
 - Pre-cut Ø25 mm, Ø30 mm, Ø35 mm and Ø40 mm: beige
- Maxi:
 - Bag capacity of +/- 710 ml (ISO 8670-2)
 - Cut-to-fit Ø15-65 mm: beige, transparent and beige with inspection window

Material:

- Soft non-woven and odor-proof plastic film
- Free of latex (PVC) and phthalates (DEHP)

Packaging:

- Boxes of 30 drainable bags and 6 filter protective stickers

Article Description	Article Number
FLEXIMA ACTIVE ROLL'UP MINI BEIGE 15-45	4620215N
FLEXIMA ACTIVE ROLL'UP MINI TRSP 15-45	4621215N
FLEXIMA ACTIVE ROLL'UP BEIGE CTF 15-60	4620115N
FLEXIMA ACTIVE ROLL'UP BEIGE 25 MM	4620125N
FLEXIMA ACTIVE ROLL'UP BEIGE 30 MM	4620130N
FLEXIMA ACTIVE ROLL'UP BEIGE 35 MM	4620135N
FLEXIMA ACTIVE ROLL'UP BEIGE 40 MM	4620140N
FLEXIMA ACTIVE ROLL'UP TRSP CTF 15-50	4621115N
FLEXIMA ACTIVE ROLL'UP MAXI BEIGE 15-65	4620380N
FLEXIMA ACTIVE ROLL'UP MAXI TRSP 15-65	4621380N