

## ACIDIC BICARBONATE CONCENTRATES AT 1+34 MIXING RATIO

Product name	K <sup>+</sup> mmol/l	Na <sup>+</sup> mmol/l	Ca <sup>++</sup> mmol/l	Mg <sup>++</sup> mmol/l	Cl <sup>-</sup> mmol/l	HCO <sub>3</sub> <sup>-</sup> mmol/l	Acetate <sup>-</sup> mmol/l	Glc. g/l	Osmolarity mOsm/l	Art. no.	Art. no.	Art. no.	Art. no.	Art. no.
										4.7 l	10 l	300 l	500 l	800 l
SW 375 A	1	138	1.25	0.5	107.5	32	3	1	289	8043	7575	7307	7230	7579
SW 376 A	1	138	1.5	0.5	108.0	32	3	1	290	8044	7576	7309	7276	39
SW 163 A	1	138	1.75	0.5	108.5	32	3	1	290	8045	3203	7310	7231	7121
SW 127 A	2	138	1.25	0.5	108.5	32	3	1	291	8046	7578	204	7278	7178
SW 380 A	2	138	1.5	0.5	109.0	32	3	1	292	8047	7580	7380	7280	7182
SW 139 A	2	138	1.75	0.5	109.5	32	3	1	292	8048	7524	7324	7224	364
SW 285 A	3	138	1.25	0.5	109.5	32	3	1	293	8049	116	96	98	3322
SW 381 A	3	138	1.5	0.5	110.0	32	3	1	294	8061	7581	7381	7281	7181
SW 178 A	3	138	1.75	0.5	110.5	32	3	1	294	8060	137	138	7234	212
SW 286 A	4	138	1.25	0.5	110.5	32	3	1	295	8059	118	93	7460	8232
SW 393 A	4	138	1.5	0.5	111.0	32	3	1	296	8057	7593	7393	7293	7298
SW 195 A	4	138	1.75	0.5	111.5	32	3	1	296	8056	131	88	8275	7123
SW 415 A	2	138	0	0.5	106.0	32	3	1	287	8071	7560			
SW 479 A	3	138	0	0.5	107.0	32	3	1	289	8072	7938			
SW 482 A	0	138	1.5	0.5	107.0	32	3	1	288	8073	7561			

Composition of the ready-to-use dialysis fluid: 1 l acidic concentrate SW xxx A and 32.775 l water for hemodialysis + 1.225 l bicarbonate concentrate 8.4% yields a ready-to-use hemodialysis solution with electrolyte concentration as mentioned above. Shelf life: 36 months

Packaging unit	Canisters p. pallet	Height (mm)	Length (mm)	Width (mm)
4.7 l canisters	128	251	188	145
10 l canisters	60	310	231	190

Packaging unit	Weight total (kg)	Height (mm)	Length (mm)	Width (mm)
300 l container	420	1,500	800	600
300 l container	424	1,440	890	605
500 l container	708.4	1,055	1,244	839
500 l container	702.2	1,115	1,244	844
500 l container	698	975	1,242	837
800 l container	1,093	1,462	1,244	840



The products mentioned may not be registered in your country or may have different dimensions or numbers of packaging units per pallet. Product information is registered in Germany. The product information as well as the art. no. may differ from country to country. Please contact your local B. Braun partner for further information.