

ACIDIC BICARBONATE CONCENTRATES AT 1+44 MIXING RATIO

Product name	K ⁺ mmol/l	Na ⁺ mmol/l	Ca ⁺⁺ mmol/l	Mg ⁺⁺ mmol/l	Cl ⁻ mmol/l	HCO ₃ ⁻ mmol/l	Acetate ⁻ mmol/l	Glc. g/l	Osmolarity mOsm/l	Art. no.	Art. no.	Art. no.	Art. no.	Art. no.
										4.7 l	10 l	300 l	500 l	800 l
SW 802 A	1	138	1.25	0.5	107.5	32	3	1	289	8026	7511	7312	7236	7124
SW 803 A	1	138	1.5	0.5	108.0	32	3	1	290	8027	7512	7354	7237	7454
SW 810 A	1	138	1.75	0.5	108.5	32	3	1	290	8031	7550	7315	7239	7125
SW 806 A	2	138	1.25	0.5	108.5	32	3	1	291	8018	7551	7316	7462	7126
SW 808 A	2	138	1.5	0.5	109.0	32	3	1	292	8019	7552	7317	7456	7455
SW 820 A	2	138	1.75	0.5	109.5	32	3	1	292	8032	7553	7318	7238	7128
SW 475 A	3	138	1.25	0.5	109.5	32	3	1	293	8021	7554	7319	7463	7129
SW 813 A	3	138	1.5	0.5	110.0	32	3	1	294	8020	7555	7320	7457	7130
SW 830 A	3	138	1.75	0.5	110.5	32	3	1	294	8033	7556	7321	7240	7131
SW 840 A	4	138	1.25	0.5	110.5	32	3	1	295	8034	7557	7322	7245	7132
SW 841 A	4	138	1.5	0.5	111.0	32	3	1	296	8035	7558	7323	7243	7133
SW 842 A	4	138	1.75	0.5	111.5	32	3	1	296	8040	7559	7325	7244	7134
SW 843 A	2	138	0	0.5	106.0	32	3	1	287	8090	7564			
SW 844 A	3	138	0	0.5	107.0	32	3	1	289	8091	7586			
SW 845 A	0	138	1.5	0.5	107.0	32	3	1	288	8092	7591			

Composition of the ready-to-use dialysis fluid: 1 l acidic concentrate SW xxx A and 42.425 l water for hemodialysis + 1.575 l bicarbonate concentrate 8.4% yields a ready-to-use hemodialysis solution with electrolyte concentrations as mentioned above. Shelf life: 36 months.

Packaging unit	Canisters p. pallet	Height (mm)	Length (mm)	Width (mm)
4.7 l canisters	128	251	188	145
10 l canisters	60	310	231	190

Packaging unit	Weight total (kg)	Height (mm)	Length (mm)	Width (mm)
300 l container	432	1,500	800	600
300 l container	436	1,440	890	605
500 l container	728.4	1,055	1,244	839
500 l container	722.2	1,115	1,244	844
500 l container	718	975	1,242	837
800 l container	1,125	1,462	1,244	840



The products mentioned may not be registered in your country or may have different dimensions or numbers of packaging units per pallet. Product information is registered in Germany. The product information as well as the art. no. may differ from country to country. Please contact your local B. Braun partner for further information.