

PRODUCT DATA SHEET

Flexima®/D` VZ^ R{ Key Flat Wafers Two-piece stoma bag systems



Flexima® Key is a flexible adhesive two-piece system with a guided coupling system for precise, simple and secure bag positioning on the wafer. The soft material of the coupling system provides flexibility which allows the wafer to adapt well to different body contours. Flexima® Key flat wafers are available in four different diameters and should be used with the Flexima® Key closed, drainable or urostomy bags.

Features

- Adhesive flexible coupling ring with guiding system allowing three different bag positions
- Available in:
 - Ø40 mm: pre-cut Ø 15 - 35 mm
 - Ø50 mm: pre-cut Ø 15 - 45 mm
 - Ø60 mm: pre-cut Ø 15 - 55 mm
 - Ø80 mm: pre-cut Ø 15 - 75 mm
- Boxes of 5 wafers

Usage

Ostomy appliance

Article Description	Article Number
SOFTIMA KEY BASE PLATE 40/15-35 MM	63040N
SOFTIMA KEY BASE PLATE 50/15-45 MM	63050N
SOFTIMA KEY BASE PLATE 60/15-55 MM	63060N
SOFTIMA KEY BASE PLATE 80/15-75 MM	63080N
SOFTIMA KEY SUPPORT PLAN PRE CUT 40/25MM	630425N
SOFTIMA KEY SUPPORT PLAN PRE CUT 50/30MM	630430N
SOFTIMA KEY SUPPORT PLAN PRE CUT 50/35MM	630436N