

PRODUCT DATA SHEET

Flexima®/D` WZ^ R® Flat Colostomy Bags One-piece stoma bag systems



With its flower shaped hydrocolloid skin protector, Flexima® standard colostomy bags adapt well to body contours and give added comfort to people with a colostomy. The colostomy bags of this range are available in three sizes, mini, midi and maxi.

Features

Colors and capacity:

- Mini:
 - Bag capacity of +/- 200 ml (ISO 8670-2)
 - Cut-to-fit Ø15-45 mm: beige
 - Pre-cut Ø20 mm, Ø25 mm, Ø30 mm and Ø35 mm: beige
- Midi:
 - Bag capacity of +/- 465 ml (ISO 8670-2)
 - Cut-to-fit Ø15-60 mm: beige and transparent
 - Pre-cut Ø25 mm, Ø30 mm, Ø35 mm, Ø40 mm and Ø45 mm: beige
- Maxi:
 - Bag capacity of +/- 590 ml (ISO 8670-2)
 - Cut-to-fit Ø15-75 mm: transparent and beige split cover

The beige split cover colostomy bags with its inspection window (maxi, cut-to-fit bags) combine both the discretion of a beige stoma bag and the visual control to check easily stoma or stool

aspect and quantity.

Highly performant degassing filter:

- The Laminar Technology (LT) filter on the mini and midi bags provides an optimized performance of degassing and deodorizing
- The High Performance (HP) filter with its high-tech breathable membrane is leak-proof and permeable to gases

Material:

- Soft non-woven and odor-proof plastic film
- Free of latex (PVC) and phthalates (DEHP)

Packaging:

- Boxes of 30 colostomy bags

Article Description	Article Number
Softima Closed Midi Beige 16-60	043016N
SOFTIMA CLOSED POUCH BEIGE 25 MM	043024N
SOFTIMA CLOSED POUCH BEIGE 30 MM	043028N
SOFTIMA CLOSED POUCH BEIGE 35 MM	043036N
SOFTIMA CLOSED POUCH BEIGE 40 MM	043040N
SOFTIMA CLOSED POUCH BEIGE 45 MM	043045N
Softima Closed Midi Transp 16-60	043616N