

AESCULAP[®] Prevision[®]

One philosophy. One platform.
More options.



AESCULAP® Prevision®

The femoral revision platform

Monobloc cemented

- 4 sizes
- 2 lengths (200 Et 240mm)
- 1 offset

Monobloc cementless

- 12 sizes
- 2 lengths (200 Et 240 mm)
- 3 offset options

Modular prox. components

- 16 sizes
- 4 diameters Et 4 lengths
- 4 offset options

Curved distal components

- 27 sizes (2 mm increments)
- 5 lengths (240 – 400 mm)
- 2 interlocking screw options

Straight distal components

- 40 sizes (1 mm increments)
- 4 lengths (200 – 320 mm)
- 2 interlocking screw options



One philosophy.

One platform.

Get more information
bbraun.com/prevision




Monobloc set configuration

- 1 basic tray (same basic tray as for modular)
- 1 additional tray



- Instruments in SQ.line® design
- Customized tray setting concept
- Intuitive & lean workflow
- Smart colour-coding for easy handling
- Intra-operative adaption of indication
- Same bed preparation for all options*



Basic modular set configuration

- 2 basic trays
- 1 additional tray** for selected option (straight or curved)
- Other indication specific optional trays

* Except for curved ** Core-sizes

AESCULAP[®] – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | www.bbraun.com