

Celsite[®] PICC-Cel

Peripherally Inserted Central Catheter



Vascular Systems

Celsite® PICC-Cel

**LATEX
FREE**

**MR
SAFE**

**300
PSI**



■ **For short to long term central venous access**
Celsite® PICC-Cel may be used for up to 3 months

■ **Drug Infusion therapy such as:**

- Antibiotherapy
- Blood Sampling and Transfusion
- Total Parenteral Nutrition (TPN)
- Antiviral Therapy
- Chemotherapy

■ **Also indicated for:**

- Central Venous Pressure Monitoring
- High Pressure Injections up to 300 psi (20.7 bar)




Grip-Lok® Securement Device for safe and clean fixation of the PICC at the skin exit



PICC-Cel Grip-Lok® can be ordered separately.
Ref. 04439010 - Box of 50 units.

Celsite® PICC-Cel Implantation Accessories – Safe and Comfortable Handling



Ordering information

Lumen	Size	Catheter Length	Guidewire Material	Guidewire Length	Reference Number
1	4 F	51 cm	Stainless Steel	Short	04434080**
				Long	04434081
			Nitinol	Short	04434082
				Long	04434083
	5 F	61 cm	Stainless Steel	Short	04434084**
				Long	04434085
2	5 F	56 cm	Stainless Steel	Short	04434088**
				Long	04434089
			Nitinol	Short	04434090
		Long		04434091	
		46 cm	Stainless Steel	Short	04434096* **
				Long	04434097*
	Nitinol		Short	04434098*	
		Long	04434099*		
	6 F	61 cm	Stainless Steel	Short	04434092**
				Long	04434093
				Nitinol	Short
			Long		04434095

* Length = 46 cm. 5 ml/sec may only be reached if the catheter is trimmed to a length of 46 cm or less. If the catheter length is 56 cm, the maximum flow rate achievable is 4 ml/sec.

** Available on request

Distributor

B. Braun Melsungen AG | Vascular Systems | Sieversufer 8 | 12359 Berlin | Germany
Phone +49 30 568207-300 | Fax +49 30 568207-130 | www.bbraun.com



B. Braun Médical | 204, avenue du Maréchal Juin | 92107 Boulogne Cedex | France
Fax +33 5 49 52 88 77

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 74 61 95-0 | Fax +49 74 61 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

The main product trademark 'Aesculap' is a registered trademark of Aesculap AG.
The product trademark „Celsite“ is a registered trademark of B. Braun Melsungen AG.

Subject to technical changes. All rights reserved.
This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.

Brochure No. 6050350

1115/0.7/2