

# PRODUCT DATA SHEET

## Flexima® Key Maxi Drainable Bags

Two-piece stoma bag systems



Flexima® Key is a flexible adhesive two-piece system with a guided coupling system for precise, simple and secure bag positioning on the wafer. The soft material of the coupling system provides flexibility which allows the wafer to adapt well to different body contours. The Flexima® Key Maxi bag presents a funnel-shaped outlet for high output.

### Features

#### Compact bag shape with large outlet:

- Ergonomic and optimized bag design for large capacity with funnel-shaped lower part
- 14 mm large outlet, easy to handle, even for people with low dexterity
- Can be connected to the Flow Collector, a 2 liter distal collecting bag
- Integrated LKP filter: leak-proof and permeable to gases
- To be used with the Flexima® Key wafer corresponding to the chosen diameter of the bag
- Boxes of 30 high-output bags

#### Easy bag positioning and removal:

- Guided positioning of the bag on the wafer thanks to the "key" system: The tab of the adhesive coupling ring on the bag is inserted into the slot of the wafer to guarantee a precise positioning
- A special cut-away section at the top of the flange facilitates the bag removal
- When the bag is removed no residue of the adhesive bag ring is left on the wafer

**Colors and capacity (ISO 8670-2):**

- Ø50 mm: capacity of +/- 780 ml, beige and transparent
- Ø60 mm: capacity of +/- 730 ml, beige and transparent
- Ø80 mm: capacity of +/- 655 ml, transparent

**Material:**

- Soft non-woven and odor-proof plastic film
- Free of latex (PVC) and phthalates (DEHP)

Article Description	Article Number
SOFTIMA KEY MAXI DRAIN POUCH BEIGE 50 MM	63550N
SOFTIMA KEY MAXI DRAIN POUCH BEIGE 60 MM	63560N
SOFTIMA KEY MAXI DRAIN POUCH TRANSP 80 M	63680N