

PRODUCT DATA SHEET

Flexima® 3S Flat Wafers

Two-piece stoma bag systems



Flexima® 3S is a 2-piece range with a guided mechanical coupling system for an easy positioning of the stoma bag onto the wafer. Flexima® 3S flat wafers are to be used with the Flexima® 3S closed, drainable or urostomy bags.

- Full skin protector layer
- Mechanical flexible coupling ring with guiding system allowing three different bag positions
- Available in diameters 45, 55, 65 and 80 mm corresponding to the bag diameters
 - Ø45 mm: cut-to-fit Ø15-30 mm and pre-cut Ø20 mm, Ø25 mm, Ø30 mm
 - Ø55 mm: cut-to-fit Ø15-40 mm and pre-cut Ø35 mm, Ø40mm
 - Ø65 mm: cut-to-fit Ø15-50 mm and pre-cut Ø45 mm
 - Ø80 mm: cut-to-fit Ø15-65 mm
 - Boxes of 5 wafers

Advantages

Flexima® 3S stand for Simple, Soft and Secure:

- **Simple** coupling system that intuitively guides you through the attachment process so that only a minimum of patient education is required. Depending on which of the three slots you choose a multiple bag orientation is possible to adapt the bag positioning to your body contours.
- **Soft** coupling system with a floating flange that requires a minimum of pressure on the abdomen to prevent discomfort during application. The coupling system remains flexible when the wafer and bag are coupled together.

- **Secure coupling system** for three reasons: Enhanced leak-proof thanks to the close fitting of the bag on the wafer, increased security through the click when wafer and bag are well locked together and finally even more fixation with the possibility to use a belt attached to one of the three hooks on the bag.

Article Description	Article Number
FLEXIMA 3S BP FLAT 45/ 15-30MM	936415N
FLEXIMA 3S BP FLAT 45/ 20MM	936420N
FLEXIMA 3S BP FLAT 45/ 25MM	936425N
FLEXIMA 3S BP FLAT 45/ 30MM	936430N
FLEXIMA 3S BP FLAT 55/ 15-40MM	936515N
FLEXIMA 3S BP FLAT 55/ 35MM	936535N
FLEXIMA 3S BP FLAT 65/ 15-50MM	936615N
FLEXIMA 3S BP FLAT 65/ 45MM	936645N
FLEXIMA 3S BP FLAT 80/ 15-65MM	936815N