

PRODUCT DATA SHEET

Flexima® 3S Urostomy Bags

Two-piece stoma bag systems



Flexima® 3S is a 2-piece range with a guided mechanical coupling system for an easy positioning of the stoma bag onto the wafer. The Flexima® 3S urostomy bags offer an anti-reflux valve to avoid urine backflow and an easy to handle and secure B. Braun drainage tap with 'open' and 'close' positions.

Secure outlet:

- Uro-specific outlet with open/closed position
- Fits in perfectly and is strengthened by a click and a final locking system
- Easy to handle and secure: no residual drops on fingers
- Can be connected to a night drainage bag by using the universal connector provided in the box (click at connection)
- The universal connector has an anti-twist design

Colors and capacity (volume as per ISO 8670-2):

- Ø45 mm: capacity of +/- 675 ml, beige with inspection window and transparent
- Ø55 mm: capacity of +/- 635 ml, beige with inspection window and transparent
- Ø65 mm: capacity of +/- 620 ml, beige with inspection window and transparent

Material and packaging:

- Soft non-woven and odor-proof plastic film
- Free of latex (PVC) and phthalates (DEHP)
- To be used with the Flexima® 3S wafer corresponding to the chosen diameter of the bag

B. Braun Melsungen AG | 34212 Melsungen | Germany
Tel. +49 5661 71-0 | Fax +49 5661 71-4567 | www.bbraun.com

- Box of 30 urostomy bags and 30 universal connectors

Advantages

Flexima® 3S stand for Simple, Soft and Secure:

- **Simple** coupling system that intuitively guides you through the attachment process so that only a minimum of patient education is required. Depending on which of the three slots you choose a multiple bag orientation is possible to adapt the bag positioning to your body contours.
- **Soft** coupling system with a floating flange that requires a minimum of pressure on the abdomen to prevent discomfort during application. The coupling system remains flexible when the wafer and bag are coupled together.
- **Secure** coupling system for three reasons: Enhanced leak-proof thanks to the close fitting of the bag on the wafer, increased security through the click when wafer and bag are well locked together and finally even more fixation with the possibility to use a belt attached to one of the three hooks on the bag.

Article Description	Article Number
FLEXIMA 3S URO BEIGE WITH WINDOW 45MM	934045N
FLEXIMA 3S URO BEIGE WITH WINDOW 55MM	934055N
FLEXIMA 3S URO BEIGE WITH WINDOW 65MM	934065N
FLEXIMA 3S URO TRP 45MM	934145N
FLEXIMA 3S URO TRP 55MM	934155N
FLEXIMA 3S URO TRP 65MM	934165N