

PRODUCT DATA SHEET

Flexima®/D` WZ^ R{ Uro Silk Flat Urostomy Bags

One-piece urostomy bag



With its flower shaped hydrocolloid skin protector, Flexima® Uro Silk bags adapt well to body contours and give added comfort to people with a urostomy. The Flexima® Uro Silk convex urostomy bags are designed for people with flush or slightly retracted urostoma. They offer an anti-reflux valve to avoid urine backflow and an easy to handle and secure B. Braun drainage tap with 'open' and 'close' positions.

Features

Secure outlet:

- Uro-specific outlet with open/closed position
- Fits in perfectly and is strengthened by a click and a final locking system
- Easy to handle, secure and hygienic: no residual drops on fingers
- Can be connected to a night drainage bag by using the universal connector provided in the box (click at connection)
- The universal connector has an anti-twist design and can freely rotate to avoid tubing twisting

Colors and bag capacity:

- Bag capacity (ISO 8670-2): +/- 630 ml
- Cut-to-fit Ø12-55 mm: transparent and beige
- Pre-cut Ø20 mm, Ø25 mm, Ø30 mm and Ø35 mm: transparent

Material and Packaging:

- Soft non-woven and odor-proof plastic film
- Free of latex (PVC) and phthalates (DEHP)
- Boxes of 30 urostomy bags and 30 universal connectors

Article Description	Article Number
SOFTIMA URO SILK TRANSP 12-55MM	043913N
SOFTIMA URO SILK BEIGE 12-55MM	043914N