



Manuplast® | Technical Data Sheet

Single Use PE Gloves (non sterile)

Device variants		Manuplast PE Gloves	
REF		6075010	6075053
Size		S (lady)	L (gentleman)
Dimension	length	285 mm ± 10 %	300 mm ± 10 %
	width total	227 mm ± 5 %	245 mm ± 5 %
	width at palm	153 mm ± 5 %	155 mm ± 5 %
	wall thickness	0.023 mm ± 10 %	
Material	low density polyethylene (LDPE) powder-free, latex-free, free of plastizisers		
Design	shape	flat, ambidextrous, welded seam	
	surface	embossed	
	colour	clear transparent	
Logistics	packaging	dispenser pack with 100 pieces	
		transportation carton with 100 packages	
	storage	room temperature, protect from humidity and sunlight	
	shelf life	5 years	
Sterilization	non sterile		
Quality certificate	ISO 9001		
Product classification	Consumer product		